

Spray Probio Tab San - 300 ml - empty

GENERAL

High-quality, reusable spray bottle for use with Probio Tab San.

- Atomizes 1ml per trigger action
- Ergonomic vaporisation thanks to Continuous Mist Spray Technology



USAGE

1. Fill the Probio Tab San Spray with 300 ml water (see marking on the label)
2. Add a dose of 1 Probio Tab San in the spray container, then seal the bottle with the spray nozzle
3. Wait for about 3 minutes until the tablet has completely dissolved
4. The spray is now ready for use

Only use this product for its specified purposes. If in doubt, first test the product in an inconspicuous area. Greenspeed is not responsible for any damages if the product is used incorrectly. Consult the Safety Data Sheet (SDS) for additional information.



1 TAB + 300 ML AQUA

Spray Probio Tab San - 300 ml - empty

TECHNICAL SPECIFICATIONS

Dimensions product (LxHxW) 8 x 25 x 5 cm

LOGISTIC DATA

Product

Article	4003642
Unit	Bottle
Dimensions	8 x 25 x 5 cm
Gross weight	95 g
Net weight	100 g
EAN code	EAN 5420015516191

Packing

Quantity per packing	110
Dimensions	60 x 46 x 40 cm
Gross weight	10500 g
EAN code	EAN 5420015516467

Pallet

Quantity per pallet	1320
Layers	3
Units per layer	440

